



SKT NEW SERIES STEREO MICROSCOPES

Quality — Crafted in Japan

Dependable — All metal construction

Priced Right — Affordable dependable
quality





SKT-2BT

FEATURES

- * Crisp, erect images over a widefield of view.
- * Binocular body is inclined at 45° for great ease of observation.
- * Binocular body can be raised or lowered on pillar stand.
- * Interpupillary distance is adjustable between 51mm and 75mm.
- * Dioptic adjustment is possible with one eyepiece tube.
- * Focusing is by rack and pinion with safety stopper at the lower end of-traverse.
- * Equipped with pairs of objectives TURRET MOUNT for rapid magnification change.
- * Magnification range is from 5X to 80X with optional eyepieces.
- * Stages are equipped with stage clips.
- * Illumination employs 12V, 10W bulb.
- * Incident and transmitted lights are designed to operate separately.
- * Supplied in styrofoam case.

Weight & Size

Model	Net wt	Size (W x D x H)
SKT 1, 2, 3	1.7 Kg	115 x 170 x 255mm
SKT 1B, 2B, 3B	2.1 Kg	115 x 170 x 260mm
SKT 1BT, 2BT, 3BT	2.5 Kg	117 x 190 x 270mm

Ordering Information (Standard specifications)

Model	Objective	Eye-piece	Base	Stage	Illuminator
SKT-1	1X & 2X	WF10X	Plain	Black/white	None
SKT-1B	1X & 2X	WF10X	w/trans.	Black/white	Incident
SKT-1BT	1X & 2X	WF10X	w/trans.	Frosted	Incident & Transmitted
SKT-2	1X & 3X	WF10X	Plain	Black/white	None
SKT-2B	1X & 3X	WF10X	w/trans.	Black/white	Incident
SKT-2BT	1X & 3X	WF10X	w/trans.	Frosted	Incident & Transmitted
SKT-3	2X & 4X	WF10X	Plain	Black/white	None
SKT-3B	2X & 4X	WF10X	w/trans.	Black/white	Incident
SKT-3BT	2X & 4X	WF10X	w/trans.	Frosted	Incident & Transmitted

Please specify mains voltage when ordering.

Optical Data

Objectives	Eyepieces	WF5X (Option)		WF10X Standard
		Magnification	5X, 10X	10X, 20X
1X & 2X	Field of View	22mm, 11mm	20mm, 10mm	
1X & 3X	Magnification	5X, 15X	10X, 30X	
	Field of View	22mm, 7.3mm	20mm, 6.5mm	
2X & 4X	Magnification	10X, 20X	20X, 40X	
	Field of View	11mm, 5.5mm	10mm, 5mm	

Objectives	Eyepieces	WF15X (Option)	WF20X (Option)	Working Distance
		15X, 30X	20X, 40X	
1X & 2X	Field of View	15mm, 7.5mm	10mm, 6.5mm	73mm
1X & 3X	15X, 45X	20X, 60X		55mm
	Field of View	15mm, 5mm	10mm, 4.3mm	
2X & 4X	30X, 45X	40X, 80X		45mm
	Field of View	7.5mm, 3.7mm	6.5mm, 3.2mm	

ACCESSORIES FOR SKT SERIES STEREO MICROSCOPES

Interchangeable Eyepieces through the all SK Series Stereo Microscopes.



MA101
WF5X, paired



MA102
WF10X, paired
(Standard)



MA103
WF15X paired




MA104
WF20X, paired



MA105
WF10X with
reticule 0.1mm


EYEPIECE MICROMETER/OBJECTIVE MICROMETER




MA523
Cross-line




MA506
10mm divided
into 100 units




MA509
5mm divided
into 100 units




MA510
10mm square
subdivided into
100 squares.



MA285 1mm divided into 100 parts



MA287 10mm square subdivided into 100 squares.



MA264



MA263

Illuminators

MA264 Versatile Illuminator with prefocused 3-lens condensing system, detachable universal arm, base transformer, 6.5V, 2.75A bulb.

MA263 Free Stand Illuminator High intensity, long life daylight, miniature flood bulb. Protected in a vented, adjustable lamp.

Please specify mains voltage when ordering.

MA264/05 Spare bulb for MA264

MA263/05 Spare bulb 120V 30W for MA263

MA263/10 Spare bulb 220/240V 30W for MA263

Other accessories

MA111 Spare bulb 12V, 10W for SKT

MA126 Wooden cabinet for SKT

Agent:



MEIJI TECHNO CO.,LTD.

6, Oi-670, Oi-machi, Iruma-gun
Saitama 356, Japan
Phone:492-67-0911 Fax:492-69-0691
Telex:25573 MEIJTCHN
Cable:MEIJITCHNO OI

MEIJI TECHNO AMERICA

2186 Bering Drive
San Jose, CA., 95131, USA
Phone:408-428-9654 Fax:408-428-0472
Toll free:800-832-0060